

LaMotte UV/VIS Spectrophotometer

UV/VIS Spectrophotometer

Order Code 7-2000-UV

- Programmable wide range spectrophotometer can accept up to 200 user curves in memory
- Pre-programmed to analyze selected tests on a direct-reading digital display
- Advanced key pad with 22 buttons simplifies data entry
- Includes 1 quartz 10mm square cuvette, 6 round 25mm 10mL test vials, universal sample holder, and square sample holder

Options:

K3 Analyst Software [Code 7-2000-UV-CD]

Wavelength range:	190-1100nm
Wavelength Accuracy:	±2nm
Wavelength Repeatability:	±1nm
Bandwidth:	4nm
Optical System:	Single Beam, grating 1200 lines/mm
Photometric Accuracy:	±0.004 A at 0.5 A
Communication Ports:	USB / RS-232
Power Requirements:	100-240V
Instrument Weight & Size:	14kg [31 lbs] / 73cm x 60cm x 34cm [19"x15"x8.5"]
Shipping Weight & Size:	17.5kg [39 lbs] / 104cm x 82cm x 60cm [26"x21"x15"]

80+ Test Factors including:

Aluminum
Ammonia
Cadmium
Chromium
Cobalt
COD LR
COD HR
Copper
Cyanide
Fluoride
Iron
Lead
Manganese LR
Manganese HR
Nickel
Nitrate
Nitrite
Phenol
Phosphate
Sulfate
Zinc



Water Testing Leader Since 1919!